



# ...for use in an Emergency

The Escape-Chair® is an indispensable solution for transporting people down the stairs in an emergency!

Example:

- Fire or smoke alarm
- Power failure
- Bomb alert
- Lift failure

## Advantages

- Ready for use in seconds
- Does not block the stairways
- You decide how quickly to descend
- Makes the evacuation flow
- Fits into every evacuation plan



### DUST COVER

Heat-resistant dust cover for storing the Escape-Chair®.



### SIGN

Luminous sign DIN 67510



### PADDED SEAT

People who are unable to sit in an upright position can be transported with the padded seat fixed onto the Escape-Chair®



### ESCAPE STRETCHER®

Flexible stretcher that helps transfer people, who may be lying on the ground or in a bed, into the Escape-Chair®.



a trademark of van Leeuwen Companies

### Life saving evacuation products

P.O. Box 2693, NL-6401 DD Heerlen  
 T +44 (0)207 2450695  
 F +44 (0)207 2450085  
 E info@escape-mobility.com  
 W www.escape-mobility.com

# The Chair







# The Chair that saves lives

Everyone uses the lifts daily, but do not realise that these lifts cannot be used if they breakdown or may not be used when there is an emergency (such as fire). Then they have to leave the premises via the emergency staircase. In practice, you find that not everyone is able to do this. Even if you have a good evacuation plan in place, it is possible that there are people present that could be at an unexpected potential risk during an emergency.

## ESCAPE CHAIR

MOBILITY COMPANY

With the help of the Escape-Chair® just one operator can bring a physically-dependent person down the stairs to safety. This evacuation chair is ready in seconds, easy to use and descends quickly and controllably down the stairs. The Escape-Chair® comes in various designs, so that for every business and situation, a solution can be found.

## ESCAPE TRAINING

MOBILITY COMPANY

We provide training in the use of our evacuation aids.

## ESCAPE ANALYSIS

MOBILITY COMPANY

We provide a building/evacuation analysis, at your premises. Here, we can provide you with advice on suitable escape routes and recommend evacuation aids to suit your situation.

## ESCAPE MAINTENANCE


MOBILITY COMPANY

So that you are at all times confident that your Escape-Chair® remains in good working order, we will send you an easy to follow checklist, so that you can service the Escape-Chair® yourself. If the glide belts need to be replaced, we will provide you with your first set of interchangeable glide belt cassettes, free of charge.

## COMFORT

### Technical information


- Weight: 13,35 kg
- Load: 136 kg
- Measurements (folded):
  - Height 120 cm
  - Width 57 cm
  - Depth 21 cm
- Two-length adjustable extending handle
- Anti-slip foam grip on handle
- Double ergonomic, curved back support
- Interchangeable glide belt cassettes
- Padded headrest
- Arm rests
- Foot rest



## STANDARDPLUS

### Technical information

- Weight: 11,2 kg
- Load: 136 kg
- Measurements (folded):
  - Height 120 cm
  - Width 51 cm
  - Depth 21 cm
- Two-length adjustable extending handle
- Anti-slip foam grip on handle
- Double ergonomic, curved back support
- Interchangeable glide belt cassettes
- Padded headrest



## STANDARD


### Technical information

- Weight: 10 kg
- Load: 136 kg
- Measurements (folded):
  - Height 102 cm
  - Width 51 cm
  - Depth 18 cm
- Standard extending handle
- Singular back support
- Interchangeable glide belt cassettes



## Carry model

In the majority of cases, the emergency escape route goes downwards. However, there are some circumstances where the escape route goes onto a higher level. For these situations it may be necessary to have an evacuation chair that can also be used as a carry chair. The Escape-Chair® Carry is fitted with carrying handles located on the top and the underside of the chair. So, with the help of two or more people, the chair can be used in this way, too.



### HOW DOES THE CHAIR BRAKE?

The unique glide/brake system makes the braking of the chair possible on different kinds of stairs, from wood to concrete. This braking system is optimised by the weight of the evacuee, together with the pressure of the operator pushing down on the handle in the direction of the stairs.

## 6 YEAR WARRANTY

Because we use long-lasting and high quality materials, the Escape-Chair® requires minimum servicing and as a result, we can guarantee a 6 year warranty.

We reserve the right to modify the technical and optical specifications.

